

DOSING GUIDE # 1

HOW

CLINICAL-GRADE

LIGHT THERAPY IS DELIVERED

And Why Guesswork Is Obsolete

LUXE

The World's First Smart
Full-Body Light Therapy Panel

luxxeredlight.com

INTRODUCTION

Red light therapy has been used in research environments for decades.

But most consumer panels were built without solving the most important variable:

Dose control.

For years, users have relied on:

- Distance estimates
- Time suggestions
- Generic positioning advice

Without real-time measurement, exposure is assumed – not verified.

This guide explains:

- What dose actually means
- Why distance changes everything
- How irradiance determines biological response
- Why visible dosing changes the future of light therapy

THE THREE VARIABLES THAT DETERMINE DOSE

In photobiomodulation, dose is not brightness.

Dose is determined by:

1. **Irradiance** (mW/cm²)
2. **Time** (minutes)
3. **Distance** (which changes irradiance)

$$\text{Dose} = \text{Irradiance} \times \text{Time}$$

If irradiance drops because you move farther away, your total delivered dose changes.

Without measuring irradiance at the body, true dosing cannot be confirmed.

THE PROBLEM WITH STATIC PANELS

Traditional panels:

- Emit a fixed output
- Assume a specific distance
- Require the user to remain still
- Provide no feedback on received irradiance

As distance changes:

Irradiance decreases exponentially.

Without adjustment, dosing becomes **variable**.

Variable dosing leads to inconsistent exposure.

In research environments, irradiance is measured.

In most home panels, it is estimated.

INTRODUCING INTELLIGENT DOSING

Luxe introduced IntelliDose™ — a patent-pending sensor system that:

- Detects user distance
- Measures received irradiance
- Adjusts LED output automatically (Mode-1 & Mode-2)
- Displays real-time dosing metrics

Instead of **assuming exposure**,
the system maintains it.

Full-body light therapy becomes **measurable**.

THE FIRST PANEL THAT SHOWS THE DOSE

Luxxe's front LCD displays:

- Distance from panel
- Irradiance reaching the body
- Total joules delivered
- Joules per square centimeter (J/cm²)
- Wavelengths in use
- Remaining session duration

This transparency allows users to:

- Understand their exposure
- Replicate clinical parameters
- Eliminate guesswork

Light therapy moves from **invisible to measurable.**

MODE-1: AUTOMATED CLINICAL PROTOCOLS

Luxxe presets are modeled after five peer-reviewed clinical studies:

- **Cardiovascular** — SPIE Proc. 11221
- **Cognitive** — PMID 37018063
- **Skin Health** — PMID 24286286
- **Pain Relief** — PMID 36359198
- **Exercise Recovery** — PMID 37099210

Mode-1 maintains target irradiance and time automatically.

No manual adjustment required.

MODE-2: DESIGNED FOR REAL-WORLD BIOLOGY

Melanin influences light absorption.

Mode-2 was intentionally developed to ensure consistent dosing across higher melanin levels.

Light therapy becomes inclusive and consistent.

PRO MODE: RESEARCH REPLICATION

PRO Mode allows:

- Independent wavelength control
- Adjustable intensity
- Real-time irradiance **visibility**
- Dose customization

Users can replicate published **clinical parameters** with visible feedback.

The variables researchers use are now available at home.

WHY THIS MATTERS

Consistent dosing supports:

- Skin tone and texture
- Muscle recovery
- Circulatory support
- Performance and vitality

Not because the panel is stronger —
but because the exposure is controlled.

CLOSING STATEMENT

Traditional panels introduced millions to light therapy.

Luxxe builds on that foundation by adding:

- Measurement
- Automation
- Transparency

When dosing becomes **intelligent**, results become consistent.